



**PLEASE ORDER BY KIT AND PART NAME.**

**SPARE PARTS AVAILABLE:**

<u>KIT No.</u>	<u>PART NAME</u>	<u>ILLUSTRATION No.</u>
WX111	Major Service Kit	4,5,6,8,9,10,12.
WX1061	Minor Service Kit	4,5,6,9,10,12.
WX499	Needle Mount Assembly	1,3,4,5,6,7,19,20.



If this instrument becomes sluggish during use, lubricate both in front of and behind the piston with NJ Phillips Lubricant.



**njphillips.com**

*Our expertise is in your hands.*

**NJ PHILLIPS PTY LTD** ABN 36 000 082 002  
**ADDRESS:** LOCKED BAG 8, GOSFORD, NSW 2250 AUSTRALIA  
**TELEPHONE:** +61 2 4340 2044 **FAX:** +61 2 4340 1991  
**EMAIL:** [njp1@njphillips.com.au](mailto:njp1@njphillips.com.au) **MADE IN AUSTRALIA**

AS OUR POLICY IS ONE OF CONTINUOUS IMPROVEMENT THE MANUFACTURER RESERVES THE RIGHT TO ALTER THESE SPECIFICATIONS AT ANY TIME. ALL PRODUCTS PRODUCED BY NJ PHILLIPS PTY LIMITED, ARE IDENTIFIED BY A UNIQUE BATCH NUMBER. THIS IDENTIFICATION NUMBER IS AFFIXED TO THE PRODUCT TO ALLOW TRACEABILITY BY THE MANUFACTURER AND MUST NOT BE REMOVED IF PRODUCT INTEGRITY IS TO BE MAINTAINED.

PAS1198 | QL686-R8